


Silicone Penrose Drain Tubes

- Penrose Drain Tube
- Penrose Straight Drain Tube
- Special Penrose Drain Tube



FORTUNE MEDICAL INSTRUMENT CORP.

6Fl., No. 29, Sec 2, Jungjeng E. Rd. Danshuei Jen, Taipei Hsien, TAIWAN
TEL: (886)-2-2624-2233 FAX: (886)-2-2624-2266
Email: fortune1@fortunemed.com info@fortunemed.com

www.FortuneMed.com

ISO 9001 & CE Certified

FDA 510(k) CE 0123 

Fortune® Silicone Penrose Drainage Tubes are used for post-operative drainage. Soft and smooth tube minimizes trauma during insertion and withdrawal, the tube draws the fluids to the exterior using the capillary phenomenon. Transparent tube with X-ray opaque line allows for easy visual inspection, fluid observation, and also confirmation of intubated tubes using X-ray; white tube has X-ray opaque effectiveness throughout the entire length. Three types of Penrose are available for selection:

- Penrose Drain Tube
- Penrose Straight Drain Tube
- Special Penrose Drain Tube

Silicone Penrose Drain Tube:

Cat. No.	O.D.	Description
2010-0405	5 mm	300 mm X-ray opaque line
2010-0406	6 mm	
2010-0408	8 mm	
2010-0410	10 mm	
2010-0412	12 mm	
2010-0415	15 mm	
2010-0419	19 mm	
2010-0425	25 mm	500 mm X-ray opaque line
2010-1405	5 mm	
2010-1406	6 mm	
2010-1408	8 mm	
2010-1410	10 mm	
2010-1412	12 mm	
2010-1415	15 mm	
2010-1419	19 mm	
2010-1425	25 mm	



Silicone Penrose Straight Drain Tube:

Cat. No.	Size	Width	Length	Description
2011-0020	S	20 mm	400 mm	X-ray opaque effectiveness
2011-0025	M	25 mm		
2011-0035	L	35 mm		
2011-0045	XL	45 mm		



Silicone Special Penrose Drain Tube:

Cat. No	Penrose Drain			Round Tube	
	O.D.	Length	X-ray	I.D. x O.D.	Length
2012-0006	6 mm	300 mm	Yes	3 x 5 mm	600 mm
2012-0008	8 mm			3 x 5 mm	
2012-0010	10 mm			4 x 6 mm	
2012-0012	12 mm			4 x 6 mm	